



Certificate of Analysis

BITUMEN MC 3000

No	Characteristic	Value	Test Method
1	Kinematic Viscosity at 60 °C (cSt)	3000-6000	ASTM D-2170
2	Flash Point (top open cup) (°C)	66 min	ASTM D-3143
3	Distillate, %Volume of Total Distillate to at 360°C: To 225°C To 260°C To 316°C	- 0-15 15-75	ASTM D-402
4	Residue from Distillation to 360°C (%Vol by difference)	80 min	ASTM D-402
Test on Residue from Distillation			
1	Viscosity at 60°C	30-120	ASTM D-2170
2	Penetration at 25°C (dmm)	120-250	ASTM D-5
3	Ductility at 25°C (cm)	100 min	ASTM D-113
4	Solubility in Trichloroethylene (%wt)	99 min	ASTM D-2042
5	Water Content (%vol)	0.2 max	ASTM D-95