



## Certificate of Analysis

### BITUMEN MC 70

No	Characteristic	Value	Test Method
1	Kinematic Viscosity@ 60°C (cSt)	70-140	ASTM D-2170
2	Flash Point (Top open cup) (°C)	38 min	ASTM D-3143
3	Distillate, %volume of total distillate to @ 360°C: to 225°C to 260°C to 316°C	20 max 20-60 65-90	ASTM D-402 ASTM D-402
	Residue from Distillation to 360°C (%Vol by difference)	55 min	
<b>Test on Residue from Distillation</b>			
1	Viscosity@ 60°C (cSt)	30-120	ASTM D-2170
2	Penetration @ 25°C (dmm)	120-250	ASTM D-5
3	Ductility @ 25°C (cm)	100 min	ASTM D-113
4	Solubility in Trichloroethylene (%wt)	99 min	ASTM D-2042
5	Water Content (%vol)	0.2 max	ASTM D-95