



Certificate of Analysis

BITUMEN MC 800

No	Characteristic	Value	Test Method
1	Kinematic Viscosity@ 60°C (cSt)	800-1600	ASTM D-2170
2	Flash Point (top open cup) (°C)	66 min	ASTM D-3143
3	Water Content (%vol)	0.2 max	ASTM D-95
4	Distillate, %volume of total distillate to @ 360°C: to 225°C to 260°C to 316°C	- 0-35 45-80	ASTM D-402
5	Residue from Distillation to 360°C (%Vol by difference)	75 min	ASTM D-402
Test on Residue from Distillation			
1	Penetration @ 25°C (dmm)	120-250	ASTM D-5
2	Ductility @ 25°C (cm)	100 min	ASTM D-113
3	Solubility in Trichloroethylene (%wt)	99 min	ASTM D-2042